

Information Overload

Managing Data in WebTMA with the Mass Import Tool





Dan Maddux

Technical Operations Specialist dan.maddux@tmasystems.com

Based in Tulsa, Dan has a 15-year career at TMA Systems, progressing from the support center to a pivotal role in Professional Services since 2010. He excels in SAML support and has a notable expertise in Microsoft Excel in addition to leading efforts in external application implementation.





Agenda

- 1. Accessing the Mass Import tool
- 2. Excel Templates
- 3. Importing Information
- 4. Reviewing/Editing information
- 5. Importing Your Data

What is Mass Import and how does it work?

Paste 🞸	Font	Alignmen	, %	ber	Conditional Format as T Cell Styles	Formatting v able v v	Cells	Editing	Sensit	tivity Add		alyze Cr Data and	reate PDF d share link
1 🗘	× ✓	fx				0							
A	в	C	D	E	F	G	н	1	,	K	L	Μ	N



Current Import Methods



Manual Input

--Proper insert/update statement IF EXISTS (SELECT 1 FROM F_UDF WHERE UDF_ITEM_FK = '1995') UPDATE F_UDF SET UDF_30CHAR2 = '010308' WHERE UDF_ITEM_FK = '1995' AND UDF_ITEMTYPE_FK = 7 ELSE INSERT INTO F_UDF (UDF_ITEMTYPE_FK, UDF_ITEM_FK, UDF_30CHAR2) VALUES ('7', '1995', '010308') --Insert Linked Documents UPDATE f_linkeddocs set lds_document = BulkColumn FROM f_linkedDocs, Openrowset (BULK 'D:\Share\DeLaSalle HighSchool \1000\1000_3_1311_area_Locker_Rm_Handicap_shower_parts.pdf', Single_Blob) Doc Where LDS_Filename = '1000_3_1311_area_Locker_Rm_Handicap_shower_parts.pdf' --backup database through query

TMA **Services**

declare @dbname nvarchar (50) = db name() declare @path nvarchar (200) = concat('''D:\Databases\backups\',@dbname,'_',(cast (Getdate()as date)),'.bak''') Professional (GetCate()as date)), .bak) declare @sql nvarchar (300) = concat('BACKUP DATABASE ',@dbname, ' TO DISK = ,@path) print @sql exec sp_sqlexec @p1 = @sql

	General	I Information			
		Code Transaction Type Connection Type	Area Import SM - Area Import Comma Delimited Text File	Active Description	✓ Area Import
		Server Name Database Name Port # User Name Password		Table Name Key Column Name Processed Column Name Error Column Name File Path	C:\Ulimport\Areas
	Mapping	gs Grid			
		TMA Column		Remote Column	
	۲	Building Code		Position: 1	
Universal	۲	Area Type Code		Position: 2	
Intorfaco	۲	Area Description		Position: 3	
Intenace	۲	Area #		Position: 4	
(IIT / IIDT)	۲	Location ID		Position: 5	
		Active?		Position: 6	
					~ /

"objectState": "Unchanged", "id": 0, "code": "string",
"alternateAccount": "string", "name": "string", "active": true, "isHeader": true, "accountTypeId": 0,
"accountTypeCode": "string", "accountTypeDescription": "string", "creatorId": 0, "createdDate": "2023-09-14T13:41:02.026Z", "modifierId": 0, "modifiedDate": "2023-09-14T13:41:02.026Z", "accountSubTypeId": 0, "accountSubTypeCode": "string", "accountSubTypeDescription": "string", "accountGroupId": 0, "accountGroupCode": "string", "accountGroupDescription": "string", "departmentId": 0, "departmentCode": "string", "departmentName": "string", "currencyId": 0, "currencyExchangeId": 0, **API Integration** "currencyExchangeDate": "2023-09-14T13:41:02.026Z", "encumbered": 0,

Mass Import Tool



© 2025 TMA SYSTEMS, LLC

Mass Import Tool

	🗰 File 🗸 1	Transactions 🛩 Material 🗸	Accounting 🗸 Organization 🗸 Ad	imin 🗸 Reports 🗸 🛧 🗸
Import				× Close
Configure Import	2 Review and Import		3 Import Complete	
DATA TYPE Select a Data Type V				
Biomed Utility Meter				
Device				
Next >			Universal Interface	
			Setup WeChat Template	
			BIM Template Setup	
			Archive and Purge	
			Mass Import >	Dashboard
				Create Import
2025 TMA SYSTEMS LLC				

Import Summary

Import summary allows you to see what has been completed so far. It is strictly a summary page.

Import Summary	
2 Imports Completed Partial	O Failed
Import List	
Biomed-25	Completed Jul 24, 2023, 6:22:46 PM
Created: 0 Updated: 1 Errored: 1	L TMA User
Biomed-24	Completed Jul 24, 2023, 6:21:30 PM
Created: 1 Updated: 1 Errored: 0	TMA User

Excel Templates

Import			+ Create Import
HISTORYImport SummaryImport List	Excel Templates		
TEMPLATES Excel Templates 	Biomed Equipment Biomed Equipment excel import template.	 Modified: Jun 19, 2023 BiomedImportTemplate.xlsx 	±

Eventually, all import templates will be listed. Following the template is key to a successful data import.



Working in Excel



Working in Excel

	A	В	С	D	E	F	G	Н		J
1	Control Id	Alternate Tag # 💌	Make Name	Biomed Type 💽 🔽	Device Type 💌	Parent Tag 🛱 💌	Facility 🔽	Building Name 📘 🔽	Floor Code 💌	Area # 🔽 🛛
2	8675		FLUKE	Defibrillators	1000		Premier City Hospital	Childrens Building		
З	309		FLUKE	Defibrillators	1000		Premier City Hospital	Childrens Building		
4										

Ensure the items you put in the spreadsheet connect to actual items in WebTMA:

- Lookups Types
- Facility
- Building
- Technician

WebTMA will not automatically add these items because they are in the spreadsheet. It will only add the biomed equipment, or whichever type of data you're importing.

Working in Excel

D	E	F	G
Biomed Type 💽 🔽	Data Reference Column 💌	Device Type 💌	Parent Tag 🛱 💌
Defibrillators	=VLOOKUP(Biomed!D2,'Data Reference	1000	
Defibrillators		1000	

Extra sheets and extra columns will need to be removed and cleaned up.



Save your file in an accessible location



Once your file is cleaned up and ready, save it in an accessible location. It is recommended that you save the file for reference in the future. However, after the import is complete, it may be discarded once data validation is finished.

Importing your data

On the Mass Import page, you will see the +Create Import.

Click there to start the import process.



Select your data type

	🗰 File 🗸 Transactions 🗸 Ma	terial 🗸 Accounting 🗸 Organization 🗸	Admin 🗸 Reports 🗸 ★	× ()
Import			× c	lose
1 Configure Import	2 Review and Imp	ort	3 Import Complete	
DATA TYPE Select a Data Type				
+ Choose			DATA TYPE	
			Select a Da	ata Type 🗸
Next >			Biomed	
			Utility Me	ter
			Device	

Types of Imports





- Region
- District
- Facility
- Building
- Floor
- Area





Asset Management

- Asset
- Biomed Equipment
- Device
- Equipment
- Tools
- Vehicles

Finance

- Accounts
- Departments



Inventory Management

- Parts
- Warehouses
- Part Inventory
- Vendors



Meters and Utilities

- Fuel and Oil
- Meters
- Meter Readings
- Meter Type
- Utility Meters
- Utility Tickets

Choose your file







Validating your file

		1 Configure Imp	ort		2 Review and I	mport		Import C	omplete	
	2 Tota	l Records		2 Records to Cre	eate	0 Records	s to Update	0 R	ecords to Skip	
< Back	Imp	port >								
									– १ 🖍	
Ū .	7 co	NTROL ID	ALTERNATE TAG #	маке наме 🍸	BIOMED TYPE	DEVICE TYPE	PARENT TAG #	FACILITY 7	BUILDING NAME Y	FLOOR
CRE	EATE 867	75		FLUKE	Defibrillators	1000		Premier City Hospital	Childrens Building	
CRE	EATE 309	9		FLUKE	Defibrillators	1000		Premier City Hospital	Childrens Building	
<				Chausian 1 to 0 of	O entries	1				>
				Showing 1 to 2 of	2 entries << <	1 > >>	25 🗸			

Modifying the import

On the Review and Import page, not only can you modify your import, you can add information that is not included in the spreadsheet.



					– द् 🖍
DATE PURCHASED	WARRANTY DATE $\ \ \square$	CREATES 7	STORES 7		CRITICAL LIFE S
Sep 20, 2023, 1:33:47 PM	Oct 20, 2025, 1:33:52 PM	true	true	\checkmark	false
		false	false	true	false

Import complete

		(1)		
	Config	jure Import		
2 item(s) succe	essfully queue	d.		
2 item(s) succ	essfully queue	d.		
2 item(s) succ	essfully queue	d.		
2 item(s) succo	essfully queue	d. Close	>	

After hitting Close, you will be taken back to the import home. You will see the end results on that page and you should now see your import in WebTMA 7.





Thank You / Q&A

Dan.Maddux@tmasystems.com | tmasystems.com

